

FINGERTIP Pulse Oximeter



USER'S MANUAL

3.3 Parameter setting:

Press the power button (> 0.5s), the oximeter will enter into parameter setting.

Settings	
Alm Setup *	
Alm	off
Beep	off
Restore	OK
Exit	

Figure 3.3.1

Settings	
Sounds Setup *	
SpO2 Alm Hi	100
SpO2 Alm Lo	90
PR Alm Hi	130
PR Alm Lo	50
+/-	+
Exit	

Figure 3.3.2

3.4 Operation

3.4.1 Install battery

Installing two AAA batteries into battery cassette in correct polarities and cover it.

WARNING: Do not attempt to recharge normal alkaline batteries, they may leak and may cause a fire or even explode.


3.4.2 Turn the Pulse Oximeter on/off

Put one of fingers into rubber hole of the Oximeter (it is best to put the finger thoroughly) with nail surface upward, then releasing the clamp.

Section 1 Safety

1.1 Instructions for the Safe Operation and Use of the Finger Pulse Oximeter

- Do not attempt to service the pulse oximeter. Only qualified service personnel should attempt any needed internal servicing.
- Do not use the oximeter in situations where alarms are required.
- Prolonged use or the patient's condition may require changing the sensor site periodically. Change sensor site and check skin integrity, circulatory status and correct alignment at least every 2 hours.
- SpO2 measurements may be adversely affected in the presence of high ambient light. Shield the sensor area (with a surgical towel, or direct sunlight, for example) if necessary.
- The following reasons will cause interference:
 - High-frequency electrosurgical
 - Placement of a sensor on an extremity with a blood pressure cuff arterial catheter, or intravascular line
 - The patient has hypotension severe vasoconstriction severe anemia or hypothermia
 - The patient is in cardiac arrest or is in shoe



Press power button to turn the Pulse Oximeter on. The oximeter will automatically be powered off when no finger in the device for longer than 16 seconds.

3.4.3. Read correspondent data from display screen.

3.4.4 Display Description of OLED

The display interface of "OLED" can rotate four directions with six different display modes after pressing the power button for less than 0.5s. It is shown as below:



- Fingernail polish or false fingernails may cause inaccurate SpO2 readings

1.2 Warnings

WARNING: EXPLOSION HAZARD--Do not use the oximeter in a flammable atmosphere where concentrations of flammable anesthetics or other materials may occur.

WARNING: Do not throw batteries in fire as this may cause them to explode.

WARNING: Do not attempt to recharge normal dry-cell batteries, they may leak. And may cause a fire or even explode.

WARNING: Do not use the pulse oximeter in an MRI or CT environment.

CAUTION: Keep the operating environment free of dust, vibrations, corrosive, or flammable materials, and extremes of temperature and humidity.

CAUTION: The battery must be taken out from the battery compartment if the device will not be used for a long time.

CAUTION: Do not operate the unit if it is damp or wet because of condensation or spills. Avoid using the equipment immediately after moving it from a cold environment to a warm, humid location.


CAUTION: Never use sharp or pointed objects to operate the front-panel key.

CAUTION:-The device shall only be used if the battery cover is closed.

CAUTION: The battery must be properly disposed according to local regulation after their use.

1.3 Definitions and Symbols

When plugging your finger into the Oximeter, your nail surface must be upward.



Declaration: Please use the medical alcohol to clean the rubber before each test and clean the tested finger with alcohol before and after the test. (The rubber inside of the Oximeter adopts medical rubber, which has no toxin, no harm, and brings no side effect such as allergy to the our skin).

3.5 Date Analysis

3.5.1 How to find the "date analysis".

Long-press on the power button, the oximeter will display pages 1 as shown figure 3.5.1 ~ 3.5.3, long-press page 1/3 and page 2/3 make the screen display the next page.



3.5.1



3.5.2




3.5.3


3.5.2 How to start a new analysis

Select the item of "record" page 3/3 and long-press its status from "off" to "on"; select "ok" when the display is shown as

figure 3.5.4. then return to page 3/3 choose "exit"; put the finger into oximeter rubber hole to start measurement as figure 3.4.3, The sign of "battery power" and "Rec" appear in turn during measuring, the maximum of recording time is 8 hours and the minimum time is 2 minutes.




3.5.4




3.5.5

3.5.3 How to see the analysis result


1. Take out your finger, wait for 8 seconds then the oximeter turn off. turn on the oximeter again. long-press to enter "data analysis: page 3/3", now the status of "records" is off and "summary graph" and "statistics" is ok, long

 The equipment type is BF

 Refer to user manual before application

%SpO2
Hemoglobin Saturation

PR BPM
Heart Rate(BPM)

 Low power indication

SN
Serial No.

Install two AAA batteries into battery cassette before covering its cover.

Plug one finger into rubber hole of the Oximeter (it is best to plug the finger thoroughly) before releasing the clamp with the nail upwards.

Press button on the front panel;

Don't tremble your finger when the Oximeter is working. Your body is not recommended on moving status.

Press the button on the front panel, if we want change display direction;

Read relevant datum from display screen.

If there is no signal input, oximeter can shut off automatically;

Please replace new batteries when OLED indicates the batteries are in low power.

Section 2 Introduction

2.1 Brief Device Description

The Pulse Oximeter, based on all digital technology, is intended for noninvasive spot-check measurement of functional oxygen

saturation of arterial hemoglobin (SpO2). Advanced DSP algorithm can minish the influence of motion artifact and improve measurement accuracy of low perfusion. The Oximeter can be used to measure human Hemoglobin Saturation anqj heart rate through finger. The product is suitable for family, hospital (including clinical use in internist/surgery, pediatrics, ect), Oxygen Bar, social medical organizations, physical care in sports and etc.

2.2 Intended use

This product is suitable for the hospital (including surgey, paediatrics, and clinical use), oxygen bar, sports health(using them before or aftersports, do not advising using them during the movement), and community health care, etc.

2.3 Product Features

- Lightweight for carrying and Easy-To-Use.
- Manually adjust the direction of interface.
- LED display, simultaneous display for testing value and bargraph.
- Low Battery voltage indicator.
- Automatically standby or sleep.

-press" "summary graph" shown in figure 3.5.8

2. Select "statistics" and long press to open the statistics page as shown figure 3.5.91. the meaning of "time" is the measuring time.it will show the date of max and Min about the spO2 and Bmp.

Note: summary graph and statistics can not be opened when the storage is empty

3.5.4 Setting in Figure 3.5.2

Select the SpO2 or PR Alarm limits, long-press will change the limits.


Select "+/- long-press will set the direction of changing the limits. "+" is increasing the number, "-" is decreasing the number.



3.5.6



3.5.7



3.5.8

2.4 Specifications

1. Display Type: OLED display

2. SpO₂:

Measurement range: 70%-99%

Accuracy: ±2% on the stage of 70%-99%, Unspecified (≤70%)

Resolution: ±1%

3. PR:

Measurement range: 30BPM-240BPM

Accuracy: ±1BPM or ±1%

4. Power Supply: 2 AAA 1.5V batteries

5. Dimension: 60*35*35mm

6. Environmental:

Operation Temperature	5°C-40°C
Storage Temperature	-10°C-40°C
Operation Humidity	15%-80%
Storage Humidity	10%-80%
Air Pressure	70-106kpa

Features

OLED display

- ◆ Two color OLED display, more display modes
- ◆ Low-power consumption, continuously four direction adjustable
- ◆ Low voltage indicator
- ◆ Automatically power off in 8 seconds when there is no signal
- ◆ Small in volume, light in weight, and convenient to carry

Section 4 Product classification, maintenance and solution

4.1 Classification

Management Class for Medical Devices:
Equipment

Anti-electric Shock Type: Internally powered
equipment

Anti-electric Shock Degree: Type BF
equipment

4.2 Maintenance and Preservation

- Replace the batteries timely when low voltage lamp is on.
- Clean the surface of fingertip oximeter before it is used to diagnose patients.
- Remove the batteries inside if you will not

Section 3 Installation, Setup, and Operation

3.1. Description of the Front Panel (as figure 3.1.1)




Figure 3.1.1 Parts of front&back panel

Table 3.1.1 Part Definition and Description

Item	Name	Descriptio
1	Power button	Turn on the machine, direction change and Parameter setting
2	OLED Panel	Display the SP02/PR data&Plethysmogram

3.2 Display

After switch on, the OLED display of the Oximeter is as follow

operate the Oximeter for a long time.


- It would be better to preserve the product in -10-40 C (14-104 °F) and humidity is 10%-80%.
- It is recommended that the product should be kept dry anytime. A wet ambience might affect its lifetime and even damage the product.
- Please follow the law of the local government to deal with used batteries.

4.3 Product declaration

Guidance and manufacture's declaration-electromagnetic radiation-
for other EQUIPMENTS and SYSTEMS

The Pulse Oximeter is designed to be used in specified electromagnetic environment . Users of the Pulse Oximeter must use it in the following environments.		
Radiation Test	Compliance	Electromagnetic environment—guidance
RF interference CISPR 11	Group 1	RF signal of Pulse Oximeter is simply created by its internal function. Therefore, its RF interference is very low and is not likely to cause any interference to nearby electronic equipment.
RF interference CISPR 11	Class B	The Pulse Oximeter applies to all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

4.4 Possible problems and effective solutions



3.3 Parameter setting:

Press the power button (> 0.5s), the oximeter will enter into parameter setting.

There are two submenu for choice:

When the '*' signal is shown on the 'Sounds setup', press the button (>0.5) and enter into the sounds setting menu (figure 3.3.1), you can press the button in turn to select the item, and press the button to change the data you need. Select or to increase or decrease the number of settings.

When the '*', signal is shown on the 'Aim Setup', press the button (> 0.5s) and enter into the alarm setting menu (Figure 3.3.2). Press the button to set the on/off for the alarm and beep.

Problem	Possible reason	Solution
SpO ₂ or PR can not be shown normally	<ol style="list-style-type: none"> 1. Finger is not plugged correctly 2. Patient's Oxyhemoglobin value is too low to be measured 	<ol style="list-style-type: none"> 1. Retry by plugging the finger 2. Try more times. If you can make sure there is no problem in the product, please go to hospital timely for exact diagnosis
SpO ₂ or PR is shown unsteady	<ol style="list-style-type: none"> 1. The finger might not be plugged deep enough 2. Finger is trembling or the patient is on movement status 	<ol style="list-style-type: none"> 1. Retry by plugging the finger 2. Please remain at rest
The Oximeter can not be turned on	<ol style="list-style-type: none"> 1. Inadequate power or power off 2. Batteries might be installed incorrectly 3. The Oximeter might be damaged 	<ol style="list-style-type: none"> 1. Please replace the batteries 2. Please reinstall the batteries 3. Please contact with local customer service centre
Indication lamps are suddenly off	<ol style="list-style-type: none"> 1. The product automatically shuts off when no signal is detected in 8 seconds 2. Inadequate power 	<ol style="list-style-type: none"> 1. Normal 2. Replace the batteries