

DemeMASK™

# Small and Childrens Face Masks

**ASTM Level 3**



Öä d ä ~ ö ä Æ ~ K



Microscope Experts since 1979

100 Lauman Lane, Suite A, Hicksville, NY 11801  
Tel: (877) 877-7274 | Fax: (516) 801-2046  
Email: [Info@nyscopes.com](mailto:Info@nyscopes.com)  
[www.microscopeinternational.com](http://www.microscopeinternational.com)



**MADE IN USA**

DemeMASK™

**Black Disposable Face Mask Small**  
**ASTM Level 3**



## Black disposable Face Mask Small

The DemeMASK disposable face mask is a fantastic protective barrier against external pathogens. It is recommended for use in common environments such as school, work or the supermarket. If the contaminants are more dangerous to health, a respirator should be used.

- DemeTECH's Black Face mask is manufactured using the highest quality of American made materials, and 100% manufactured in The United States.
- The DemeTECH Black Mask is made with the same materials as other DemeTECH Surgical Mask.
- Adjustable to the nose with a mouldable aluminum strap and rubber grip for greater comfort, without squeezing or leaving marks on the skin. All four edges of the mask are heat sealed for superior tear resistance.



Black Disposable Face Mask Small

Product Number	Type	Valve	Nose Clip	Strap Style	Strap Material	Box Qty
DT-MSK-003B-SMALL	Disposable	NO	YES	Earloops	Latex Free Elastic	25 Units



### Technical specifications

- ASTM Level 3
- 3 protective layers filter >98% Particulates
- The edges of the mask are heat sealed
- 2 adjustment bands
- The high-density polypropylene outer layer provides strong resistance against liquids
- The soft polypropylene inner layer provides comfort for long periods of use



DemeMASK™

**Pink Disposable Face Mask Small**  
**ASTM Level 3**



## Pink disposable Face Mask Small

The DemeMASK disposable face mask is a fantastic protective barrier against external pathogens. It is recommended for use in common environments such as school, work or the supermarket. If the contaminants are more dangerous to health, a respirator should be used.

- DemeTECH's Pink Face mask is manufactured using the highest quality of American made materials, and 100% manufactured in The United States.
- The DemeTECH Pink Mask is made with the same materials as other DemeTECH Surgical Mask.
- Adjustable to the nose with a mouldable aluminum strap and rubber grip for greater comfort, without squeezing or leaving marks on the skin. All four edges of the mask are heat sealed for superior tear resistance.



Pink Disposable Face Mask Small

Product Number	Type	Valve	Nose Clip	Strap Style	Strap Material	Box Qty
DT-MSK-PNS	Disposable	NO	YES	Earloops	Latex Free Elastic	50 Units



### Technical specifications

- ASTM Level 3
- 3 protective layers filter >98% Particulates
- The edges of the mask are heat sealed
- 2 adjustment bands
- The high-density polypropylene outer layer provides strong resistance against liquids
- The soft polypropylene inner layer provides comfort for long periods of use

DemeMASK™

**Green Disposable Face Mask Small**  
**ASTM Level 3**



## Green disposable Face Mask Small

The DemeMASK disposable face mask is a fantastic protective barrier against external pathogens. It is recommended for use in common environments such as school, work or the supermarket. If the contaminants are more dangerous to health, a respirator should be used.

- DemeTECH's Green Face mask is manufactured using the highest quality of American made materials, and 100% manufactured in The United States.
- The DemeTECH Green Mask is made with the same materials as other DemeTECH Surgical Mask.
- Adjustable to the nose with a mouldable aluminum strap and rubber grip for greater comfort, without squeezing or leaving marks on the skin. All four edges of the mask are heat sealed for superior tear resistance.



Green Disposable Face Mask Small

Product Number	Type	Valve	Nose Clip	Strap Style	Strap Material	Box Qty
DT-MSK-GRS	Disposable	NO	YES	Earloops	Latex Free Elastic	50 Units



### Technical specifications

- ASTM Level 3
- 3 protective layers filter >98% Particulates
- The edges of the mask are heat sealed
- 2 adjustment bands
- The high-density polypropylene outer layer provides strong resistance against liquids
- The soft polypropylene inner layer provides comfort for long periods of use

DemeMASK™

**Purple Disposable Face Mask Small**  
**ASTM Level 3**



## Purple disposable Face Mask Small

The DemeMASK disposable face mask is a fantastic protective barrier against external pathogens. It is recommended for use in common environments such as school, work or the supermarket. If the contaminants are more dangerous to health, a respirator should be used.

- DemeTECH's Purple Face mask is manufactured using the highest quality of American made materials, and 100% manufactured in The United States.
- The DemeTECH Purple Mask is made with the same materials as other DemeTECH Surgical Mask.
- Adjustable to the nose with a mouldable aluminum strap and rubber grip for greater comfort, without squeezing or leaving marks on the skin. All four edges of the mask are heat sealed for superior tear resistance.



Purple Disposable Face Mask Small

Product Number	Type	Valve	Nose Clip	Strap Style	Strap Material	Box Qty
DT-MSK-PRS	Disposable	NO	YES	Earloops	Latex Free Elastic	50 Units



### Technical specifications

- ASTM Level 3
- 3 protective layers filter >98% Particulates
- The edges of the mask are heat sealed
- 2 adjustment bands
- The high-density polypropylene outer layer provides strong resistance against liquids
- The soft polypropylene inner layer provides comfort for long periods of use

DemeMASK™

**Yellow Disposable Face Mask Small**  
**ASTM Level 3**

## Yellow disposable Face Mask Small

The DemeMASK disposable face mask is a fantastic protective barrier against external pathogens. It is recommended for use in common environments such as school, work or the supermarket. If the contaminants are more dangerous to health, a respirator should be used.

- DemeTECH's Yellow Face mask is manufactured using the highest quality of American made materials, and 100% manufactured in The United States.
- The DemeTECH Yellow Mask is made with the same materials as other DemeTECH Surgical Mask.
- Adjustable to the nose with a mouldable aluminum strap and rubber grip for greater comfort, without squeezing or leaving marks on the skin. All four edges of the mask are heat sealed for superior tear resistance.



Yellow Disposable Face Mask Small

Product Number	Type	Valve	Nose Clip	Strap Style	Strap Material	Box Qty
DT-MSK-YLS	Disposable	NO	YES	Earloops	Latex Free Elastic	50 Units



### Technical specifications

- ASTM Level 3
- 3 protective layers filter >98% Particulates
- The edges of the mask are heat sealed
- 2 adjustment bands
- The high-density polypropylene outer layer provides strong resistance against liquids
- The soft polypropylene inner layer provides comfort for long periods of use