

DemeMASK™

Surgical MASK

FDA Clearance: K201479



DemeMASK™

Surgical Mask MSK-003

FDA Clearance: K201479

Level 3 Surgical Mask

ASTM Level 3 Masks feature a highly protective three layer construction, giving you protection, breathability, and comfort in the OR.

THREE LAYER FILTER PROTECTION



**Non-Woven
Protective Layer**

Blue outer layer ensures proper wear, and provides high protection against fluids, high breathability, and excellent protection against airborne bacteria.

Filter Layer

Provides at least 98% filtration against airborne particles and bacteria.

Inner Comfort Layer

Enables wearers to comfortably wear mask for extended periods of time.


High Protection across all devices

Providing Your Staff With Superb Comfort, Comfort, without sacrificing reliable protection

At DemeTECH, we continually work with excellent quality materials, every mask is manufactured with the wearers protection in mind. The result is a highly protective surgical mask, with excellent comfort and world class protection.

Comfort

- Superb breathability – While filtering harmful bacteria, our masks allow for easy air flow in both directions, do the wearer stays cool and comfortable.



Bi-Directional Air Flow

Innovation


- Unique Device coding system to ensure staff knows which mask to use



A revolutionary new way to fight viruses

Protection

- Filtration Efficiency: DemeTECH masks meet the bacterial filtration efficiency and differential pressure standards set forth in "ASTM F2100 Standard Specification for the performance of Materials Used in Medical Facial Masks"
- Fluid Protection - Choose from three levels of fluid resistance: 80mmHg, 120mmHg, and 160mmHg.



High Levels of fluid resistance against blood and other airborne bodily fluids

Öä dā` çā/å` K

 **New York
Microscope
Company**
Microscope Experts since 1979

100 Lauman Lane, Suite A, Hicksville, NY 11801
Tel: (877) 877-7274 | Fax: (516) 801-2046
Email: Info@nyscopes.com
www.microscopeinternational.com