

DISPOSABLE ISOLATION GOWN





Isolation Gowns are used to protect patients from possible infection that may be sourced in hospitals. Non-surgical gowns are Class I devices intended to protect the wearer from the transfer of microorganisms and body fluids in low or minimal risk patient isolation situations.

1	COMPOSITON	% 100 PP SB Non Woven Fabric TNT
2	FABRIC WEIGHT	30-40-50 gsm
3	COLOR	White-Blue-Other
4	SIZE EU	S, M,L,XL,XXL
5	MADE IN	Türkiye
6	CARE INSTRUCTION	Single Use



DIMENSIONS in cm



	S	M	L	XL	XXL
A	145,0	150,0	155,0	160,0	160,0
B	110,0	115,0	120,0	123,0	130,0
C	19,0	20,0	21,0	22,0	23,5
D	56,0	56,0	57,0	58,0	64,0
E	54,0	54,0	56,0	56,0	60,0
F	32,0	32,0	31,0	32,0	36,0
G	60,0	60,0	60,0	60,0	60,0
H	5,0	5,0	5,0	5,0	5,0

PACKING

Isolation Gowns are packed in PE bags. Packages show the information below if requested:

- Product series name & product name
- Set type (dimensions and type of drape)
- Set contents with quantities
- Manufacturing date
- Lot number
- Reference number
- Warnings (don't use open packages, for single use only, latex free)
- Manufacturer contact & company information
- CE Label

QUALITY CERTIFICATE & COMPLIANCE

Isolation Gowns are CE Class I products. All process across Farabi is run in compliance with EN ISO 9001 & EN 14126 CE certificate is available.

FDA registration No: **10068316**



SAFETY INSTRUCTIONS

- Non Hazardous
- Store in cool and dry area
- No skin irritation
- In case of rashes discontinue use wash hands with plenty of water.