

20000 SERIES

POWDER FREE, CLEAR VINYL EXAMINATION GLOVES

Our most economical exam glove that provides durable protection for routine medical tasks. Our glove is DEHP free and latex free, soft, comfortable, and is an effective barrier to biological pathogens.



GLOVE SPECIFICATIONS

WEIGHT & DIMENSION			
Reorder#	Size	Length	Palm Width
20210	X-Small	250mm	75 ± 5mm
20211	Small	250mm	85 ± 5mm
20212	Medium	250mm	95 ± 5mm
20213	Large	250mm	105 ± 5mm
20214	X-Large	250mm	115 ± 5mm

MIL THICKNESS (SINGLE WALL)		
	Palm	Finger
All Sizes	3.8 ± 0.2 mil	4.2 ± 0.2 mil

PHYSICAL PROPERTIES	
	Before Aging
Tensile Strength Mpa (minimum)	13
Ultimate Elongation % (minimum)	380

PERFORMANCE REQUIREMENTS			
Characteristic	Related Defect	Inspection Level	AQL
Freedom from holes	Holes, lumps	G - 1/ S - 4 Multiple Normal	1.5
Physical dimension	Length, width and thickness	S - 2	1.5
Physical properties	Before aging and after accelerated aging	S - 2	1.5
Visual defect	Stains etc.	S - 4	1.5

PACKAGING	
10 boxes/case	100 gloves/box (by weight)

New York Microscope Company
Phone: (877) 877-7274
Email: info@nyscopes.com
www.microscopeinternational.com

⚠ WARNING: This product contains cancer and reproductive harm chemicals. Visit www.P65Warnings.ca.gov/product.