

10000 SERIES

POWDER FREE LATEX EXAMINATION GLOVES – TEXTURED

Our most popular powder free exam glove provides exceptional comfort and quality at a competitive price point without sacrificing consistent dependability. This glove has a micro-textured grip and a precise fit providing reduced hand fatigue.



GLOVE SPECIFICATIONS

WEIGHT & DIMENSION			
Reorder#	Size	Length	Palm Width
10110	X-Small	238 - 242mm	75 ± 5mm
10111	Small	238 - 242mm	85 ± 5mm
10112	Medium	238 - 242mm	95 ± 5mm
10113	Large	238 - 242mm	105 ± 5mm
10114	X-Large	238 - 242mm	115 ± 5mm

MIL THICKNESS (SINGLE WALL)		
	Palm	Finger
All Sizes	4.0 mil	4.5 mil

PHYSICAL PROPERTIES	
	Before Aging
Tensile Strength Mpa (minimum)	18
Ultimate Elongation % (minimum)	650

PERFORMANCE REQUIREMENTS			
Characteristic	Related Defect	Inspection Level	AQL
Freedom from holes	Holes, lumps	G - 1 / S - 4 Multiple Normal	1.5
Physical dimension	Length, width and thickness	S - 2	4.0
Physical properties	Before aging and after accelerated aging	S - 2	4.0
Visual defect	Stains etc.	S - 4	2.5

PACKAGING	
10 boxes/case	100 gloves/box (by weight)