

16000 SERIES

POWDER FREE LATEX EXAMINATION GLOVES – TEXTURED - 12” EXTENDED CUFF

Our premium heavy gauge latex exam glove with a 12 inch long cuff. Designed for use in high risk situations by First Responders, Fire Departments, Police, Ambulance EMTs, Coroners, and Morticians. This glove is available in 15 mil, 12 mil, and 10 mil thicknesses, allowing high risk professionals to select the proper thickness to meet their safety needs. This glove is blue color, fully textured, with an extended cuff for maximum protection.



GLOVE SPECIFICATIONS

WEIGHT & DIMENSION			
Reorder#	Size	Length	Palm Width
16011	Small	300 ± 10mm	80 ± 10mm
16012	Medium	300 ± 10mm	95 ± 10mm
16013	Large	300 ± 10mm	110 ± 10mm
16014	X-Large	300 ± 10mm	≥ 110mm
16015	XX-Large	300 ± 10mm	≥ 113mm

MIL THICKNESS (SINGLE WALL)		
	Palm	Finger
All Sizes	13 mil	15 mil

PHYSICAL PROPERTIES	
	Before Aging
Tensile Strength Mpa (minimum)	18
Ultimate Elongation % (minimum)	650

PERFORMANCE REQUIREMENTS			
Characteristic	Related Defect	Inspection Level	AQL
Free from holes	Holes	G - 1	1.5
Visible defects	Stains, lumps, holes, weight	S - 4	2.5
Dimensions	Width, length, thickness, weight	S - 2	4.0
Physical properties	Before & after aging	S - 2	4.0

PACKAGING	
10 boxes/case	50 gloves/box (by weight)