



30000 SERIES

POWDER FREE NITRILE EXAMINATION GLOVES FULLY TEXTURED – BLUE

Our premium, 100% LATEX FREE, powder free Nitrile Exam Glove provides excellent chemical splash resistance and puncture resistance. This glove provides exceptional elasticity and strength, is comfortable to wear for extended use and has a textured grip for use in wet and dry conditions.



GLOVE SPECIFICATIONS

WEIGHT & DIMENSION			
Reorder#	Size	Length	Palm Width
30310	X-Small	238 – 242mm	75 ± 5mm
30311	Small	238 – 242mm	85 ± 5mm
30312	Medium	238 – 242mm	95 ± 5mm
30313	Large	238 – 242mm	105 ± 5mm
30314	X-Large	238 – 242mm	115 ± 5mm

MIL THICKNESS (SINGLE WALL)		
	Palm	Finger
All Sizes	4.5 mil	5.0 mil

PHYSICAL PROPERTIES	
	Before Aging
Tensile Strength Mpa (minimum)	14
Ultimate Elongation % (minimum)	500

PERFORMANCE REQUIREMENTS			
Characteristic	Related Defect	Inspection Level	AQL
Freedom from holes	Holes, lumps	G - 1/ S - 4 Multiple Normal	1.5
Physical dimension	Length, width and thickness	S - 2	4.0
Physical properties	Before aging and after accelerated aging	S - 2	4.0
Visual defect	Stains etc.	S - 4	4.0

PACKAGING	
10 boxes/case	100 gloves/box (by weight)



100 Lauman Lane, Suite A, Hicksville, NY 11801
Tel: (877) 877-7274 | Fax: (516) 801-2046
Email: info@nyscopes.com
www.microscopeinternational.com