

30350/100 SERIES

LIGHT WEIGHT POWDER FREE NITRILE EXAMINATION GLOVES

Our light weight, powder free, nitrile examination gloves provide superior comfort, dexterity and protection. These gloves are fully micro-textured for better grip. They offer enhanced tactile sensitivity and provide biological barrier and offers splash resistance to most chemical compounds.



GLOVE SPECIFICATIONS

WEIGHT & DIMENSION

Reorder#	Size	Length	Palm Width
30351/100	Small	240 ± 2mm	85 ± 5mm
30352/100	Medium	240 ± 2mm	95 ± 5mm
30353/100	Large	240 ± 2mm	110 ± 5mm
30354/100	X-Large	240 ± 2mm	115 ± 5mm

MIL THICKNESS (SINGLE WALL)

	Palm	Finger
All Sizes	3.0 mil	3.5 mil

PHYSICAL PROPERTIES

	Before Aging
Tensile Strength Mpa (minimum)	14
Ultimate Elongation % (minimum)	500

PERFORMANCE REQUIREMENTS

Characteristic	Related Defect	Inspection Level	AQL
Freedom from holes	Holes, lumps	FDA As per 21 CFR Vol. 8 / G-1 Single Normal	1.5
Physical dimension	Length, width and thickness	S - 2	4.0
Physical properties	Before aging and after accelerated aging	S - 2	4.0
Visual defect	Stains etc.	S - 4	2.5

PACKAGING

10 boxes/case	100 gloves/box, XL-95gloves/box (by weight)
---------------	---