



18000 SERIES

POLY SILK – POWDER FREE, BLACK LATEX EXAMINATION GLOVES
– FULLY TEXTURED



The premium black latex exam glove – providing the soft, supple, powder free characteristics for maximum comfort and extended wear. The gloves are fully textured with superior grip for dry or wet applications.



GLOVE SPECIFICATIONS

WEIGHT & DIMENSION			
Reorder#	Size	Length	Palm Width
18111	Small	240 ± 2mm	85 ± 5mm
18112	Medium	240 ± 2mm	95 ± 5mm
18113	Large	240 ± 2mm	105 ± 5mm
18114	X-Large	240 ± 2mm	115 ± 5mm

MIL THICKNESS (SINGLE WALL)		
	Palm	Finger
All Sizes	3.9 mil	4.3 mil

PHYSICAL PROPERTIES	
	Before Aging
Tensile Strength Mpa (minimum)	18
Ultimate Elongation % (minimum)	650

PERFORMANCE REQUIREMENTS			
Characteristic	Related Defect	Inspection Level	AQL
Freedom from holes & lumps	Holes, lumps	G - 1 / S - 4 Multiple Normal	1.5
Physical dimension	Length, width and thickness	S - 2	4.0
Physical properties	Before aging and after accelerated aging	S - 2	4.0
Visual defect	Stains etc.	S - 4	2.5

PACKAGING	
10 boxes/case	100 gloves/box (by weight)

New York Microscope Company
Phone: (877) 877-7274
Email: info@nyscopes.com
www.microscopeinternational.com